


Product Datasheet

Her2 (ERBB2) antibody (orb1322802)

Description	Mouse monoclonal antibody to Her2 (ERBB2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 676-1255 of human ERBB2(NP_004439) produced in HEK293T cell
Target	ERBB2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	137.7 kDa
Uniprot ID	P04626
NCBI	NM_004448
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...