


Product Datasheet

PVRL1 (NECTIN1) antibody (orb1322807)

Description	Mouse monoclonal antibody to PVRL1 (NECTIN1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PVRL1(NP_002846) produced in HEK293T cell.
Target	NECTIN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54 kDa
Uniprot ID	Q15223
NCBI	NM_002855
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...