


Product Datasheet

DTYMK antibody (orb1322812)

Description	Mouse monoclonal antibody to DTYMK
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human DTYMK(NP_036277) produced in HEK293T cell.
Target	DTYMK
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23.6 kDa
Uniprot ID	P23919
NCBI	NM_012145
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...