


Product Datasheet

CTBP1-DT antibody (orb1322823)

Description	Mouse monoclonal antibody to CTBP1-DT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-186 of human C4orf42(NP_443093) produced in E.coil.
Target	C4orf42
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	19.5 kDa
NCBI	NM_052861
Dilution Range	IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent
staining of COS7
cells...