

Product Datasheet

Cytochrome P450 17A1 (CYP17A1) antibody (orb1322850)

Description Mouse monoclonal antibody to Cytochrome P450 17A1 (CYP17A1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human CYP17A1(NP_000093) produced in HEK293T cell.

Target CYP17A1

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 57.2 kDa

Uniprot ID [P05093](#)

NCBI [NM_000102](#)

Dilution Range WB 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...