


Product Datasheet

CCM2 antibody (orb1322862)

Description	Mouse monoclonal antibody to CCM2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCM2(NP_113631) produced in HEK293T cell.
Target	CCM2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	48.7 kDa
Uniprot ID	Q9BSQ5
NCBI	NM_031443
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...