


Product Datasheet

IFI35 antibody (orb1322907)

Description	Mouse monoclonal antibody to IFI35
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human IFI35(NP_005524) produced in HEK293T cell.
Target	IFI35
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	31.6 kDa
Uniprot ID	P80217
NCBI	NM_005533
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...