


Product Datasheet

Nudel (NDEL1) antibody (orb1322941)

Description	Mouse monoclonal antibody to Nudel (NDEL1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human NDEL1(NP_001020750) produced in HEK293 cell.
Target	NDEL1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	36.9 kDa
Uniprot ID	Q9GZM8
NCBI	NM_001025579
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...