


Product Datasheet

HSPBP1 antibody (orb1322942)

Description	Mouse monoclonal antibody to HSPBP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human HSPBP1(NP_036399) produced in HEK293T cell.
Target	HSPBP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	39.1 kDa
Uniprot ID	Q9NZL4
NCBI	NM_012267
Dilution Range	WB 1:500, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...