


Product Datasheet

TRAP alpha (SSR1) antibody (orb1322976)

Description	Mouse monoclonal antibody to TRAP alpha (SSR1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human SSR1 (NP_003135) produced in HEK293T cell.
Target	SSR1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.1 kDa
Uniprot ID	P43307
NCBI	NM_003144
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...