


Product Datasheet

OBFC1 (STN1) antibody (orb1322995)

Description	Mouse monoclonal antibody to OBFC1 (STN1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human OBFC1 (NP_079204) produced in HEK293T cell.
Target	STN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	41.9 kDa
Uniprot ID	Q9H668
NCBI	NM_024928
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...