


Product Datasheet

PSMD2 antibody (orb1323025)

Description	Mouse monoclonal antibody to PSMD2
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PSMD2(NP_002799) produced in HEK293T cell.
Target	PSMD2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	100 kDa
Uniprot ID	Q13200
NCBI	NM_002808
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...