


Product Datasheet

Tropomodulin 1 (TMOD1) antibody (orb1323044)

Description	Mouse monoclonal antibody to Tropomodulin 1 (TMOD1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TMOD1(NP_003266) produced in HEK293T cell.
Target	TMOD1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	40.4 kDa
Uniprot ID	P28289
NCBI	NM_003275
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...