


Product Datasheet

CNDP2 antibody (orb1323049)

Description	Mouse monoclonal antibody to CNDP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human CNDP2(NP_060705) produced in HEK293T cell.
Target	CNDP2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.7 kDa
Uniprot ID	Q96KP4
NCBI	NM_018235
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...