


Product Datasheet

GDAP1L1 antibody (orb1323059)

Description	Mouse monoclonal antibody to GDAP1L1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human GDAP1L1 (NP_076939) produced in HEK293T cell.
Target	GDAP1L1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	41.8 kDa
Uniprot ID	Q96MZ0
NCBI	NM_024034
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...