


Product Datasheet

LMCD1 antibody (orb1323071)

Description	Mouse monoclonal antibody to LMCD1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human LMCD1(NP_055398) produced in HEK293T cell.
Target	LMCD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	40.7 kDa
Uniprot ID	Q9NZU5
NCBI	NM_014583
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...