


Product Datasheet

LSM1 antibody (orb1323099)

Description	Mouse monoclonal antibody to LSM1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-133 of human LSM1 (NP_055277) produced in E.coli.
Target	LSM1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	15 kDa
Uniprot ID	O15116
NCBI	NM_014462
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...