

Product Datasheet

Spermine synthase (SMS) antibody (orb1323109)

Description Mouse monoclonal antibody to Spermine synthase (SMS)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human SMS (NP_004586) produced in HEK293T cell.

Target SMS

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 41.1 kDa

Uniprot ID [P52788](#)

NCBI [NM_004595](#)

Dilution Range WB 1:2000, IF 1:50~100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...