


Product Datasheet

SM22 alpha (TAGLN) antibody (orb1323112)

Description	Mouse monoclonal antibody to SM22 alpha (TAGLN)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human TAGLN(NP_003177) produced in HEK293T cell.
Target	TAGLN
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22.4 kDa
Uniprot ID	Q01995
NCBI	NM_003186
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...