


Product Datasheet

ILVBL antibody (orb1323123)

Description	Mouse monoclonal antibody to ILVBL
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ILVBL(NP_006835) produced in HEK293T cell.
Target	ILVBL
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	67.7 kDa
Uniprot ID	A1L0T0
NCBI	NM_006844
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...