


Product Datasheet

XPF (ERCC4) antibody (orb1323136)

Description	Mouse monoclonal antibody to XPF (ERCC4)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human ERCC4 (NP_005227) produced in HEK293T cell.
Target	ERCC4
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	104.3 kDa
Uniprot ID	Q92889
NCBI	NM_005236
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...