


Product Datasheet

beta IV Tubulin (TUBB4A) antibody (orb1323194)

Description	Mouse monoclonal antibody to beta IV Tubulin (TUBB4A)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TUBB4 (NP_006078) produced in HEK293T cell.
Target	TUBB4A
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	49.4 kDa
Uniprot ID	P04350
NCBI	NM_006087
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...