

Product Datasheet

LIM Kinase 1 (LIMK1) antibody (orb1323197)

Description Mouse monoclonal antibody to LIM Kinase 1 (LIMK1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human LIMK1(NP_002305) produced in HEK293T cell.

Target LIMK1

Note For research use only

Isotype IgG2a

Clonality Monoclonal


MW 72.4 kDa

Uniprot ID [P53667](#)

NCBI [NM_002314](#)

Dilution Range WB 1:1000, IHC 1:150, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...