


Product Datasheet

FGF21 antibody (orb1323218)

Description	Mouse monoclonal antibody to FGF21
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-209 of human FGF21(NP_061986) produced in E.coli.
Target	FGF21
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	19.4 kDa
Uniprot ID	Q9NSA1
NCBI	NM_019113
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...