


Product Datasheet

SNAP25 antibody (orb1323279)

Description	Mouse monoclonal antibody to SNAP25
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SNAP25 (NP_003072) produced in HEK293T cell.
Target	SNAP25
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23.2 kDa
Uniprot ID	P60880
NCBI	NM_003081
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...