


Product Datasheet

CDK2 antibody (orb1323342)

Description	Mouse monoclonal antibody to CDK2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human CDK2 (NP_001789) produced in HEK293T cell.
Target	CDK2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	P24941
NCBI	NM_001798
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...