


Product Datasheet

PRKAR1B antibody (orb1323355)

Description	Mouse monoclonal antibody to PRKAR1B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human PRKAR1B (NP_002726) produced in HEK293T cell.
Target	PRKAR1B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	42.9 kDa
Uniprot ID	P31321
NCBI	NM_002735
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...