


Product Datasheet

FLYWCH2 antibody (orb1323373)

Description	Mouse monoclonal antibody to FLYWCH2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human FLYWCH2 (NP_612448) produced in HEK293T cell.
Target	FLYWCH2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.4 kDa
Uniprot ID	Q96CP2
NCBI	NM_138439
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...