


Product Datasheet

PDHA1 antibody (orb1323416)

Description	Mouse monoclonal antibody to PDHA1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human PDHA1 (NP_000275) produced in HEK293T cell.
Target	PDHA1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	40.2 kDa
Uniprot ID	P08559
NCBI	NM_000284
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...