


Product Datasheet

TBC1D4 antibody (orb1323423)

Description	Mouse monoclonal antibody to TBC1D4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TBC1D4 (NP_0055647) produced in HEK293T cell.
Target	TBC1D4
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	146.4 kDa
Uniprot ID	O60343
NCBI	NM_014832
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...