


Product Datasheet

ACOT12 antibody (orb1323432)

Description	Mouse monoclonal antibody to ACOT12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human ACOT12 (NP_570123) produced in HEK293T cell.
Target	ACOT12
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	61.9 kDa
Uniprot ID	Q8WYK0
NCBI	NM_130767
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...