


Product Datasheet

SULT1C2 antibody (orb1323439)

Description	Mouse monoclonal antibody to SULT1C2
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SULT1C2(NP_001047) produced in HEK293T cell.
Target	SULT1C2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	34.7 kDa
Uniprot ID	O00338
NCBI	NM_001056
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...