


Product Datasheet

DEF8 antibody (orb1323446)

Description	Mouse monoclonal antibody to DEF8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human DEF8 (NP_997397) produced in HEK293T cell.
Target	DEF8
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	58.5 kDa
Uniprot ID	Q6ZN54
NCBI	NM_207514
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...