

Product Datasheet

Guanylate kinase (GUK1) antibody (orb1323456)

Description Mouse monoclonal antibody to Guanylate kinase (GUK1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IHC, WB

Immunogen Full length human recombinant protein of human GUK1 (NP_000849) produced in HEK293T cell.

Target GUK1

Note For research use only

Isotype IgG2b

Clonality Monoclonal


MW 21.5 kDa

Uniprot ID [Q16774](#)

NCBI [NM_000858](#)

Dilution Range WB 1:500~2000, IHC 1:150, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...