


Product Datasheet

NDUFA7 antibody (orb1323459)

Description	Mouse monoclonal antibody to NDUFA7
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human NDUFA7 (NP_004992) produced in HEK293T cell.
Target	NDUFA7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	12.4 kDa
Uniprot ID	O95182
NCBI	NM_005001
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...