


Product Datasheet

UBXN2B antibody (orb1323463)

Description	Mouse monoclonal antibody to UBXN2B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human UBXN2B (NP_001071087) produced in HEK293T cell.
Target	UBXN2B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	36.9 kDa
Uniprot ID	Q14CS0
NCBI	NM_001077619
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...