


Product Datasheet

HDHD1A (PUDP) antibody (orb1323489)

Description	Mouse monoclonal antibody to HDHD1A (PUDP)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human HDHD1A (NP_036212) produced in HEK293T cell.
Target	PUDP
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23.5 kDa
Uniprot ID	Q08623
NCBI	NM_012080
Dilution Range	WB 1:1000~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...