


Product Datasheet

BCAR1 antibody (orb1323496)

Description	Mouse monoclonal antibody to BCAR1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human BCAR1 (NP_055382) produced in HEK293T cell.
Target	BCAR1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	93.2 kDa
Uniprot ID	P56945
NCBI	NM_014567
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...