


Product Datasheet

GBA3 antibody (orb1323555)

Description	Mouse monoclonal antibody to GBA3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-150 and 370-469 of human GBA3 (NP_066024) produced in HEK293T cell.
Target	GBA3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.5 kDa
Uniprot ID	Q9H227
NCBI	NM_020973
Dilution Range	WB 1:500~2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...