


Product Datasheet

KCTD14 antibody (orb1323561)

Description	Mouse monoclonal antibody to KCTD14
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human KCTD14 (NP_076419) produced in HEK293T cell.
Target	KCTD14
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.4 kDa
Uniprot ID	Q9BQ13
NCBI	NM_023930
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...