


Product Datasheet

PKC epsilon (PRKCE) antibody (orb1323664)

Description	Mouse monoclonal antibody to PKC epsilon (PRKCE)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PRKCE (NP_005391) produced in HEK293T cell.
Target	PRKCE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	83.5 kDa
Uniprot ID	Q02156
NCBI	NM_005400
Dilution Range	WB 1:1000~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...