


Product Datasheet

HAGHL antibody (orb1323676)

Description	Mouse monoclonal antibody to HAGHL
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human HAGHL (NP_115680) produced in HEK293T cell.
Target	HAGHL
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.1 kDa
Uniprot ID	Q6PII5
NCBI	NM_032304
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...