


Product Datasheet

RTF2 antibody (orb1323701)

Description	Mouse monoclonal antibody to RTF2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human C20orf43 (NP_057491) produced in HEK293T cell.
Target	RTF2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	Q9BY42
NCBI	NM_016407
Dilution Range	WB 1:1000~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...