


Product Datasheet

CHCHD5 antibody (orb1323703)

Description	Mouse monoclonal antibody to CHCHD5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-110 of human CHCHD5(NP_115685) produced in E.coli.
Target	CHCHD5
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	12.2 kDa
Uniprot ID	Q9BSY4
NCBI	NM_032309
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...