


Product Datasheet

PIH1D2 antibody (orb1323713)

Description	Mouse monoclonal antibody to PIH1D2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PIH1D2 (NP_620144) produced in E.coli.
Target	PIH1D2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.8 kDa
Uniprot ID	Q8WWB5
NCBI	NM_138789
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...