


Product Datasheet

VCAM1 antibody (orb1323715)

Description	Mouse monoclonal antibody to VCAM1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human VCAM1(NP_001069) produced in HEK293T cell.
Target	VCAM1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	78.7 kDa
Uniprot ID	P19320
NCBI	NM_001078
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...