


Product Datasheet

KATNAL1 antibody (orb1323730)

Description	Mouse monoclonal antibody to KATNAL1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human KATNAL1 (NP_115492) produced in HEK293T cell.
Target	KATNAL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55.2 kDa
Uniprot ID	Q9BW62
NCBI	NM_032116
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...