


Product Datasheet

NUDT6 antibody (orb1323744)

Description	Mouse monoclonal antibody to NUDT6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human NUDT6 (NP_009014) produced in HEK293T cell.
Target	NUDT6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	P53370
NCBI	NM_007083
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...