


Product Datasheet

ACY3 antibody (orb1323775)

Description	Mouse monoclonal antibody to ACY3
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ACY3 (NP_542389) produced in HEK293T cell.
Target	ACY3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	35.1 kDa
Uniprot ID	Q96HD9
NCBI	NM_080658
Dilution Range	WB 1:500~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...