


Product Datasheet

Butyrylcholinesterase (BCHE) antibody (orb1323781)

Description	Mouse monoclonal antibody to Butyrylcholinesterase (BCHE)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human BCHE (NP_000046) produced in HEK293T cell.
Target	BCHE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.2 kDa
Uniprot ID	P06276
NCBI	NM_000055
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...